



Molecular Size & Mass Kit

BACKGROUND:

Students will use indirect measurements to estimate the size and mass of a molecule of oleic acid. Please note that this lab is more of a measurement exercise than a fact-finding mission. The object is to teach students about a method of measuring things they cannot see, and to give them a rough idea of the size and mass of a single molecule.

Oleic acid has the formula $C_{18}H_{34}O_2$ and the structure below:



One end of the molecule is hydrophilic (attracted to water) and the other is hydrophobic. When oleic acid is dropped on the surface of water, the molecules stand up on end, like people standing in a crowd. The molecules form a monolayer – a layer that is one molecule thick.

In this lab, oleic acid is provided in a solution of ethanol. This is done to allow only a very small volume of acid to be put on the water surface. Otherwise, you'd need a very, very large tray to do the lab.

To find the mass of a molecule, students will need to know the acid's density. If they had access to pure oleic acid, they could find the density by measuring the mass of a known volume. But since they only have access to the dilute solution, they will be provided with the density, 0.87 g/cm^3 .

Students will model the experiment first by using spherical bb's. Dispense the bb's so that each group has about 50mL of bb's.

KIT CONTENTS:

1L Oleic Acid Solution (5mL acid in Ethanol) 1 piece of chalk 12 eye droppers ~4000 bb's (~400mL)

Also needed for each group:

Large flat pan – pizza pans are ideal. Metric Ruler 100mL graduated cylinder 10mL graduated cylinder Micrometer (optional) Access to a balance

PREREQUISITES:

Students should have a good understanding of the concept of density and geometric formulas for area and volume. They should be able to solve simple algebraic formulas and solve a proportion. A good grasp of scientific notation is also very helpful.

FORMULAS AND HINTS:

Depending on the level of your students, you might want to put some or all of the following reminders on the board or discuss and practice them before the lab. High school chemistry students, for instance, should be familiar with these basic geometry formulas and algebraic manipulations. Give just enough guidance that students don't get caught in the geometry and lose track of the meaning of the lab.

- 5. (and 15) Radius is half of diameter.
- 6. Area of a circle: A = πR^2 (and 17) Cylinder volume formula, solved for height of cylinder: $h = \frac{V}{r}$
- 7. (and 18) Hint: How many bb diameters (molecule lengths) equals the height of the cylinder from #6 (#17)?
- 9. Volume of a sphere: $V = \frac{4}{3}\pi R^3$
- 10. Don't forget to subtract the mass of the container that holds the bb's.

12. (and 21) Density formula, solved for mass: $m = \rho V$

24. Percent error:

$$\% error = \frac{(actual - measured)}{actual} *100\%$$

ANSWERS:

Naturally, as in any experimental lab, you can expect a range of results. Here are the accepted values:

BB radius: 0.45cm BB mass: 0.37g Acid molecule length: 1.8 x 10⁻⁷cm Acid molecule mass: 4.68 x 10⁻²²g

Acceptable experimental results for the acid molecule will be in the correct order of magnitude (10^{-7} or 10^{-22}). If students are much farther off than that, they might want to recalculate or repeat the experiment. (Be sure to clean the pan thoroughly between trials.) Don't let students worry if they get a very high percent error (100% or more) for #23. Keep in mind that the percent difference between 1.0×10^{-7} cm and 2.0×10^{-7} cm is 100%, but they can both be considered very accurate for this lab.

REFERENCES:

Robinson, Paul. *Conceptual Physics Laboratory Manual:* 4th edition. Prentice Hall. p. 173 and p. 175

Oleic acid structure diagram: http://animalscience.tamu.edu/nutr/202s/LectureOutlines/lipids.html. (viewed 9/9/2003).

Visual acuity information:

<u>http://www.nde-ed.org/EducationResources/CommunityCollege/PenetrantTest/Introduction/visualacuity.html</u>. (viewed 9/9/2003).

ACKNOWLEDGEMENTS:

Many thanks to Mike Horton, Perris High School, for his input on this product.

RELATED PRODUCTS:

Micrometer (P1-1040). 25mm x 0.01mm.

Basic Lab Supply Pack (C3-3000). Includes beakers, graduated cylinders, flasks, test tubes, and more for outfitting student chemistry lab kits.



PO Box 2750 ANN ARBOR, MI 48106 T 800-367-6695 WWWARBORSCICOM ©2009 ARBOR SCIENTIFIC ALL RIGHTS RESERVED

USING DIRECT MEASUREMENT TO FIND MOLECULAR SIZE AND MASS

You will be indirectly measuring the size and mass of a single molecule of oleic acid. First, the process will be modeled with spherical bb's.

OBJECTIVES:

To estimate the length of a molecule of oleic acid.

To estimate the mass of a molecule of oleic acid.

MATERIALS:

BB's 100mL graduated cylinder Metric Ruler Oleic Acid solution (5mL acid in 1000mL of solution) Eye dropper Large tray or pizza pan Water Chalk dust 10mL graduated cylinder

PROCEDURE:

Model: Diameter of a BB

- 1. Use the 100mL graduated cylinder to find the volume of bb's. _____
- 2. 1mL = 1 cm³. Record the volume of bb's in cm³. _____
- 3. Pour the bb's out onto the tray and arrange them into a closely packed circle, 1 bb thick.
- 4. Measure the diameter of the circle in 3 places and find the average diameter.

Avg. Diameter: _____

5. Find the area of the circle.

Area: _____

- 6. The circle of bb's is really a very short cylinder, 1bb high. The volume of a cylinder can be found by multiplying its top area by its height. $V = A \bullet h$. Find the height of the cylinder.
- 7. What is the diameter of a single bb?
- 8. If you have a micrometer, directly measure the diameter of a bb._____

Molecular Size & Mass

Name: _____

MODEL: MASS OF A BB

9. Use the diameter from #7 to find the volume of a single bb.

Volume =_____

10. Measure the mass of the entire volume of bb's.

11. Calculate the density (ρ) of the bb's (use the total volume of bb's from #1).

$$\rho = \frac{m}{V} =$$

12. Use the density from #11 and the volume from #9 to find the mass of a single bb.

EXPERIMENT: HEIGHT OF AN ALEIC ACID MOLECULE

An oleic acid molecule is not a sphere – it is elongated like a hot dog. The acid is less dense than water. One end is attracted to water and the other is repelled, so when a drop of acid is placed on water, the molecules stay on top of the water and group side-by-side, like people in a crowd.

- 13. Fill the tray with water about 1cm deep and place it on a level surface. Allow the water to settle. Lightly scrape a piece of chalk over the water until a fine dust is on the surface. (The dust should be barely visible.)
- 14. Use the eyedropper to gently add a single drop of oleic acid solution to the center of the tray. The alcohol in the solution will dissolve into the water, but the acid will stay on top. Measure the diameter of the circle in at least 3 places and find the average.

Avg. Diameter: _____

15. Find the area of the circle.

Area: _____

PO Box 2750 ANN ARBOR MI 48106 T 800-367-6695 WWW.ARBORSCI.COM ©2009 ARBOR SCIENTIFIC ALL RIGHTS RESERVED

20. If you had access to pure oleic acid, you could measure the mass of a certain volume to find its density. Since you only have the dilute solution, use a density of 0.87 g/cm3.

19. Use the height from #19 to find the volume of a single molecule. (Approximate the molecule as a cylinder, with a diameter equal to the one-tenth of the height.)

What is the height of a single acid molecule? 18.

Experiment: Mass of a molecule of Oleic Acid

V=_____

5mL_acid _____ Volume_acid_in_1_drop 1000mL_solution Total_volume of 1 drop (#16b)

c. The acid solution contains 5mL of acid in 1000mL of solution. Use the proportion below to find the volume of acid in one drop.

The volume you found in #16c is the volume of the cylinder of acid. Find the height of the

Volume of 1 drop _____

b. Divide 3cm3 by the number of drops in 3cm3 (#16a) to find the volume of one drop.

a. Use the eye dropper to fill the 10mL graduated cylinder to a volume of 3cm3. (Remember to measure to the bottom of the meniscus.) Record the number of drops in 3cm3. Do this 3 times and find the average.

16. The volume of oleic acid in the single drop of dilute solution is much less than the volume of the whole drop. You need to know the volume of only the acid.

Molecular Size & Mass

17.

cylinder. ()

Name:

Drops in 3cm³



Name: _____

21. Use the density from #20 and the volume from #19 to find the mass of a single molecule.

Extension: Checking Accuracy

22. From 12 inches away, the human eye can see objects as small as 0.00886cm. How many oleic acid molecules lined up end-to-end would it take for a human to see the length of the line?

23. The molar mass of oleic acid is 282.4654 grams. Divide this mass by Avogadro's number (6.02×1023) to find the actual mass of one molecule.

24. Compare your experimental mass (#21) to the answer to #23 and calculate your percent error.