

## NEULOG RESISTANCE LOGGER SENSOR GUIDE



### NeuLog resistance logger sensor NUL-249

The NeuLog resistance sensor can be used for any science experiment which utilizes resistance readings. It is used in the fields of Physics, Electronics, Chemistry, Biology, etc.

The sensor comes pre-calibrated so you can start experimentation right out of the box using any of the following guides.

Hundreds of possible experimental subjects that can be performed with the NUL-249 sensor are: Electrical components resistance measurements, electrical parallel circuits, electrical series circuits and electronics.

The resistance sensor's measurement units are:

Ohm ( $\Omega$ ): The ratio between the voltage on and the resistance through an electrical component.

### Included with the sensor:

- NeuLog General Guide
- Red (positive probe) and black (negative probe) wires connected directly to the sensor's body.
- 4 mm red and black plugs for the positive and the negative wires.

### Sensor specifications

|                                  |   |
|----------------------------------|---|
| <b>Range and operation modes</b> | 0 to 200 k $\Omega$   |
| <b>ADC resolution</b>            | 15 bit  |
| <b>Resolution</b>                | (0 to 10K $\Omega$ ) 0.001 K $\Omega$<br>(10 to 60K $\Omega$ ) 0.01 K $\Omega$<br>(60 to 200K $\Omega$ ) 0.1 K $\Omega$ |
| <b>Max sample rate (S/sec)</b>   | 100   |

**Experiment Duration:** 1 second to 31 days.

### Sensor features:

- Fully digital data.
- Rugged plastic ergonomic case.
- 4 mm power and ground plugs for easy connectivity
- Push button switch for Start/Stop experiments in off line mode.
- LED indicator of experiment status (blinks while collecting data).
- Pre-calibrated sensing equipment.

**Note:** NeuLog products are intended for educational use.

## NEULOG RESISTANCE LOGGER SENSOR GUIDE



### Technical background:

The philosophy behind NeuLog's plug and play technology is based on each sensor's ability to store its own data due to an internal flash memory chip and micro-controller in each plastic NeuLog body. This technology allows the sensor to collect and then store the digital data in the correct scientific units ( $^{\circ}\text{C}$ ,  $^{\circ}\text{F}$ , Lux, %, ppm, for example).

The sensor is pre-calibrated at the factory. The built-in software in the logger can be upgraded for free at any time using the provided firmware update.

The NeuLog resistance sensor drives current to the tested resistance through a resistor in series. The sensor measures the current through the resistance and the voltage on it. Using ohm's law ( $R = V/I$ ) the resistance is calculated from the current and voltage values.

### Maintenance and storage:

- Never submerge the NeuLog plastic body in any liquid.
- Do not allow liquid into the resistance sensor's body.
- After using the probe, wipe off all excess material, liquid or residue from the sensor.
- Store in a box at room temperature out of direct sunlight.

### Warranty:

We promise to deliver our sensor free of defects in materials and workmanship. The warranty is for a period of 3 years from the date of purchase and does not cover damage of the product caused by improper use, abuse, or incorrect storage. Sensors with a shelf life such as ion selective probes have a warranty of 1 year. Should you need to act upon the warranty, please contact your distributor. Your sensor will be repaired or replaced.

**Thank you for using NeuLog!**



Flexible, simple, fast, forward thinking.

W: [www.neulog.com](http://www.neulog.com)

E: [info@neulog.com](mailto:info@neulog.com)

A: 850 St Paul Street, Suite 15, Rochester, NY 14605

P: 1.866.553.8536

V2015.5