

Pascal's Demonstrator

P1-2190



BACKGROUND:

The Pressure Syringe (aka Pascal's Demonstrator) consists of a vertical glass bulb about 3.5cm in diameter with small openings along one circumference and a branch tube containing a tight-fitting piston. It is designed to show that pressure in a fluid is transmitted equally in all directions.

OPERATION:

In short, the apparatus is filled with water, and when the piston is pushed in, a stream of water is projected from each of the openings with the same force, showing that the pressure is transmitted in all directions.

- 1. Fill a large beaker or other container with water.
- 2. Immerse the bulb in the water and, starting with the piston at its lowest position, pull the piston upward. Atmospheric pressure will force the water into the openings and fill up the bulb.
- 3. A few rapid jerks of the piston may be necessary to completely fill the bulb and the cylinder.
- 4. When the apparatus is filled it may be quickly taken from the water and the piston pressed against the liquid, forcing streams of liquid out in all directions.
- 5. The apparatus may be held in any position to produce this result, but perhaps the horizontal position is best. In this position, students will observe that the streams of water are projected about the same distance from the bulb and that, therefore, is a rough measure of the amount of force, thus verifying that the liquid pressure is transmitted in all directions.

Alternatively, the apparatus may be filled by removing the piston from the apparatus and filling from a faucet. However, with this method, some water will be lost before you can begin the demonstration. Using a vessel as described will allow you to being the demonstration immediately after filling.

RELATED PRODUCTS:

Bernoulli Bags (P6-7350) Use Bernoulli's principle to your advantage when you fill an 8-foot bag with only one breath of air!

Elasticity of Gases (P1-2075) The perfect apparatus for verifying Boyle's Law and other properties of gases. **Atmospheric Pressure Cups** (P1-2005) In this economical version of the classic Magdeburg Hemispheres, students push two rubber cups together and find that it takes tremendous force to pull them apart.